







# DeZire Mask Family



## Product Series: DeZire® Full Face (Oronasal) Masks

MODEL	NIV Mask				ePAP Face Mask	Lite CPAP Face Mask		Lite-2 Bi-level CPAP Face Mask	
									
	No Port	With Port**	No Port	With Port**	With Port*	With Port*	With Port*	With Port*	With Port* & **
REF Codes	ZV12 (S) ZV11 (M) ZV10 (L)	ZV17 (S) ZV16 (M) ZV15 (L)	ZV22 (S) ZV21 (M) ZV20 (L)	ZV27 (S) ZV26 (M) ZV25 (L)	ZE12 (S) ZE11 (M) ZE10 (L)	ZC12 (S) ZC11 (M) ZC10 (L)	ZC22 (S) ZC21 (M) ZC20 (L)	ZB22 (S) ZB21 (M) ZB20 (L)	ZB27 (S) ZB26 (M) ZB25 (L)
Cushion	Double layer (Transparent)		Single layer (Yellow)		Single layer (Yellow)	Double layer (Transparent)	Single layer (Yellow)	Single layer (Yellow)	
Anti-asphyxia valve	N/A		N/A		Yes	Yes	Yes	Yes	
Exhalation vent hole	N/A		N/A		N/A (Pressure Spring Valve)	Yes (Blue)	Yes (Blue)	No (White Cap)	
* 2 Ports for CO <sub>2</sub> / pressure monitoring (4mm)	N/A		N/A		Yes	Yes	Yes	Yes	Yes
** Port for Intubation /Feeding or Suction Tube ( elastic 3.2mm)	N/A	Yes	N/A	Yes	N/A	N/A	N/A	N/A	Yes
Headgear	Reusable		Economy		Economy	Reusable	Economy	Economy	

ABOUT DEZIRE MASK  
[www.deziremask.com](http://www.deziremask.com)  
 DeZire® is a registered trademark of GaleMed® Corporation.  
 All rights reserved.

CONTACT US  
 Email: [info@galemed.com](mailto:info@galemed.com)  
 10F, No. 109, Min-chuan E. Road, Sec. 6, Taipei, Taiwan 114  
 Tel: +886 (0) 28791-9348

WEBSITES  
 GaleMed - [www.galemed.com](http://www.galemed.com)  
 DeZire – [www.deziremask.com](http://www.deziremask.com)