


Babi.Plus Bubble CPAP System



GaleMed®

Babi Plus Bubble CPAP System



Babi Plus Bubble PAP valve

- For neonatal <10 kg
- Gas flow 1-12 liters per minute
- Continuous Positive Airway Pressure 0-10 cmH₂O

PAP Adjustment Lever
Simple pressure level changes.

Clear Water Level Marks
Allows quick set up and system maintenance.

Embedded Diffuser
Generates uniform bubbles.

Swivel Connector
In compliance with all standard 15M and 22F connectors.

Water Adjustment Port
Quick water adjustment without interruption of therapy and disconnection of expiratory circuit.

Mounting Bracket
Secures fastening of equipment on standard ISO mounts, preventing water spills.

**Bigger Bubble
Bigger Amplitude**

Babi Plus bubble CPAP 2.0 is Coming Soon

- Generate a broader frequency oscillation wave
- Re-expand the alveoli
- Maintain airway patency

