


Babi.Plus

Neonate Lung Simulator

- Anatomical design
- Better consistency in compliance
- Two different lungs simulation
- Various resistance available
- Standard ISO connections



nTest Lung



- AN0004** nTest Lung, R145x2
- AB0051** nTest Lung Module, R90
- AB0052** nTest Lung Module, R600
- AB0053** nTest Lung Module, R300
- AB0054** nTest Lung Module, R145
- AB0055** nTest Lung Module, R90

QTY: 10 pcs/box, 4 box/carton

Specifications	With Bracket	With Bracket
Tidal volume	25 ml/bag	60 ml/bag
Compliance (approx.)	1.0 ml/cmH ₂ O per bag	NA
Weight (approx.)	43 g	29 g
Dimension (W × D × H)	125×105×35 mm	
Resistance (inner diameter)	90 (3.5 mm) , 145 (3.0 mm), 300(2.5 mm) or 600 (2.0 mm) cmH ₂ O/L/s AT 1 L/s flow	
Interface connector	Y connector: 15/10 mm, 10 mm I.D., 10 mm I.D. Stranght Connector: 10 mm O.D.	
Pressure Monitoring Port	4.3 mm I.D.	

