

Babi.Plus™

nBag Neonatal Resucitator

- A safer pulmonary-assist device for premature and newborn infant up to 5kg.
- Biocompatible, latex free and Phthalate free silicone material not only reduce the risk of allergy or infections but also provide better patient feedback.
- Minimize patient risk during operation through patented design and two pressure limit levels, along with unique bag design to deliver smaller stroke volume for premature infants.
- Highly efficient and flexible oxygen reservoir system provide convenient and easy operation.



Features

Pressure Monitoring Port



Pressure Relief Valve



25 cmH₂O



40 cmH₂O

360° Swivel




20 ml



50 ml

Resuscitation Volume

Expandable Reservoir



Order Information:



REF	Descriptions	Packaging: 12 pcs/box
BR30	nBag Manual Resuscitator	
BR32	nBag Manual Resuscitator with Mask (size 0)	
Specifications: a. Adjustable 25 cmH ₂ O and 40 cmH ₂ O pressure relief valve b. Allergy-free silicone mask (size 0)		
BL10	nTest Lung Set, R145	
GB10	GiO™1 Digital Pressure Gauge	