



nMask

Lightweight and Comfort-fit Silicone Mask

The nMask is an interface to use between premature or newborn infant and:

- manual self-rebounding resuscitator
- 2 non-self-inflation resuscitator
- gas powered resuscitator
- 4 aerosol generating device
- 6 breathing circuit or anesthesia circuit



Designed Specifically for Neonates

- Provides smaller size range to reduce leakage and conform better for resuscitation and gas transfer to neonates
- Connector: 15mm O.D.
- Material: Translucent soft silicone
- Product not manufactured with material which contains DEHP and does not contain natural rubber latex
- Features round contour, ergonomic fit and rolled cushion for optimal fit
- Six cushion sizes: 25mm, 30mm, 35mm, 40mm, 45mm, 50mm
- Choose an appropriately size face mask to cover the mouth and nose for positive pressure ventilation
 Mask should not cover the eyes or overlap the chin



REF	Dscription	QTY
AN0038	nMask, 25 mm	5 pcs/box
AN0039	nMask, 30 mm	5 pcs/box
AN0040	nMask, 35 mm	5 pcs/box
AN0041	nMask, 40 mm	5 pcs/box
AN0042	nMask, 45 mm	5 pcs/box
AN0043	nMask, 50 mm	5 pcs/box